

About



Well Biotech Co., Ltd is a Sino-American joint venture dedicated to developing, manufacturing and distributing In Vitro diagnostic Assay products. These products include rapid Oral Fluid HIV1/2 antibody test for the detection of antibodies to HIV virus, Oral Fluid H Treponemal antibody test for syphilis, Oral Fluid HCV antibody test, Saliva Alcohol Test, Oral Fluid drug of abuse Screen Device, and other rapid immunoassays.

Mission



Well Biotech is committed to provide high quality diagnostic products to meet the global healthcare needs, and to empower patients to make critical decision about their health.

Quality



Well Biotech complies with the latest quality systems and employs state of art manufacturing process and cutting edge technologies. Our manufacture and research facility is TUV ISO 13485:2016 certified and China FDA authorized. Most of our products are either CFDA approved or CE registered.

Well Biotech Co., Ltd.

4F, C1 Building, No.9 Changyang Rd.,
Changzhou, Jiangsu P.R. China, 213149

Well Bioscience Inc.

9924 Mesa Rim Rd, San Diego, CA 92121

T + 86-519-8019 1660
F + 86-519-8019 1663
E sales@wellbioscience.com

www.wellbioscience.com



EC REP

Medwheat Tech Service GmbH
Max-Planck-Straße 4 85609
Aschheim b. München Germany



Orawell[®]

Rapid Diagnostics

PRODUCT CATALOG

Well Biotech Co., Ltd.



Drugs of Abuse

Product Name	Specimen	Cutoff
Oralwell® Oral Fluid Drug Screen Device	Saliva	Any Combination of the following drug assay with requested cutoff
AMP/BAR/BZO/COC/KET/MET/MDMA/MTD/MOR/OXY/PPX/THC/Alcohol/FET/TRA/6-AM/K2	Saliva	AMP25/BZO5/COC20/MET25/OPI10/THC20/MTD AMP50/BZO10/COC20/MET50/OPI40/THC40/OXY40/ng/ml 0.02% B.A.C.
6-Acetyl Morphine (6-AM)	Saliva	10 ng/ml
Cotinine (COT)	Saliva	25 or 50 ng/ml
Ketamine (KET)	Saliva	50 ng/ml
Phencyclidine (PCP)	Saliva	10 ng/ml
THC	Saliva	20 or 40 ng/ml

Drug Screen Dipcard	Urine	Single or Combination of the following drug assay with requested cutoff
6-Acetyl Morphine (6-AM)	Urine	10 ng/ml
Amphetamine (AMP)	Urine	300/500/1000 ng/ml
Barbiturates (BAR)	Urine	300 ng/ml
Benzodiazepine (BZO)	Urine	300 ng/ml
Cocaine (COC)	Urine	150/300 ng/ml
Kratom	Urine	300 ng/ml
Ethyl Glucuronide (EtG)	Urine	500/1000 ng/ml
EDDP	Urine	100/300 ng/ml
Fentanyl (FEN)	Urine	200 ng/ml
Ketamine (KET)	Urine	500/1000 ng/ml
Ecstasy (MDMA)	Urine	500 ng/ml
Methamphetamine (MET)	Urine	500/1000 ng/ml
Methadone (MTD)	Urine	300 ng/ml
Opiates (MOR)	Urine	300/1000 ng/ml
Oxycodone (OXY)	Urine	100 ng/ml
Phencyclidine (PCP)	Urine	25 ng/ml
Propoxyphene (PPX)	Urine	300 ng/ml
Tricyclic Antidepressants (TCA)	Urine	1000 ng/ml
Marijuana (THC)	Urine	50 ng/ml
Tramadol (TRA)	Urine	300 ng/ml
Synthetic Marijuana (K2)	Urine	25 ng/ml
Zolpidem (ZOL)	Urine	50/100 ng/ml
Mitragynine	Urine	1000 ng/ml



Infectious Diseases

Product Name	Specimen	Sensitivity
Oralwell® Oral Fluid HIV1/2 Antibody Screen Test	Saliva	99.2%
Oralwell® Oral Fluid HCV Antibody Screen Test	Saliva	99.2%
HIV1/2 Antibody Screen Test	WB/S/P	99.8%
HCV Antibody Screen Test	WB/S/P	99.6%
Dengue IgG & IgM Test	WB/S/P	98%
Dengue NS1 Ag Test	WB/S/P	96%
Syphilis (TP antibody) Test	WB/S/P	99.6%
Malaria	WB/S/P	50 parasites
Adenovirus	Feces	10 ⁵ vp/ml
Rotavirus	Feces	10 ⁴ cfu/ml
H. Pylori Ag/Ab	Feces	97%



Pregnancy & Fertility

Product Name	Specimen	Sensitivity
LH Ovulation Test	Urine	25 mIU/ml
FSH Infertility Test	Urine	25 mIU/ml
Pregnanediol-3 Glucuronide (PdG)	Urine	50-4000 ng/ml
Etrone-3 Glucuronide (E1G)	Urine	5-30 ng/ml
β-HCG	Urine	3-20000 mIU/ml
β-HCG	WB/S/P	5-50000 mIU/ml



Quantitative Fluorescent Immunoassay

Product Name	Specimen	Range
NT-ProBNP	WB/S/P	50-10000 pg/ml
PCT	WB/S/P	0.1-200 ng/ml
CRP	WB/S/P	0.5-200 µg/ml
SAA	WB/S/P	5-300 µg/ml
IL6	WB/S/P	3-4000 pg/ml
TSH	Serum/Plasma	0.1-5.0 µIU/ml
Troponin I	Serum/Plasma	0.1-50 ng/ml
Cortisol	Saliva	0.1-30 ng/ml



Veterinary Test

Product Name	Specimen	Sensitivity
CCV Ag Test	Feces	98%
CPV Ag Test	Feces	98%
Progesterone Test	milk	10 ng/ml



Well Biotech

Oralwell® Oral Fluid HIV 1/2 is SFDA approved rapid assay for detection of HIV 1/2 antibodies in oral fluid and delivers results in 10 minutes.

